

Features	DataViewer	Prova Free	Prova Subscription
Data analysis			
Supports fitting with: <ul style="list-style-type: none"> • Hertzian Contact • Oliver & Pharr • JKR 	✓	✓	✓
DMA analysis	✓	Improved	Improved Most filtering options
Guided data quality assessment	✓	✓ Improved	✓ Improved
Adhesive properties quantification		✓	✓
Save/load analysis projects		✓	✓
Shareable analysis SOPs			✓
Visualization			
Visualize results of full plates at once			✓
Group and compare datasets			✓
Visualization options include: <ul style="list-style-type: none"> • Heatmaps • Boxplots • Histograms • Line plots with error bars 	Limited	✓	✓ Most options
Reporting			
PDF report with traceability data			✓
Export graphs as PNG, SVG, CSV		✓	✓
Bulk CSV export	✓	✓ Improved	✓ Improved
Additional Features			
All future functionality			✓