

Prova feature list



Features	Prova Free	Prova Subscription	DataViewer
Data analysis			
Supports fitting with: <ul style="list-style-type: none"> • Hertzian Contact • Oliver & Pharr • JKR 	✓	✓	Available in DataViewer
DMA analysis	✓	✓	Improved
Guided data quality assessment	✓	✓	Improved
Quantify adhesive properties	✓	✓	New
Save/load analysis projects	✓	✓	New
Shareable analysis SOPs		✓	New
Visualization			
Visualization options include: <ul style="list-style-type: none"> • Heatmaps • Boxplots • Histograms • Line plots with error bars 	✓	✓	Improved (including new features)
Visualize results of full plates at once		✓	New
Group and compare datasets		✓	New
Reporting			
Export graphs as PNG, SVG, CSV	✓	✓	New
Bulk CSV export	✓	✓	Improved
PDF report with traceability data		✓	New
Additional Features			
All future functionality		✓	